
# gojji<sup>®</sup> Plus

Multi-Functional Monitoring System (Blood Glucose and  $\beta$ -Ketone)  
*Sistema de monitoreo multifuncional (Glucosa Sanguínea y  $\beta$ -cetona)*

## Owner's Manual

### *Manual del Propietario*

**Use of this device on multiple patients may lead to transmission of Human Immunodeficiency Virus (HIV), Hepatitis C Virus (HCV), Hepatitis B Virus (HBV), or other bloodborne pathogens.**

## Dear Gojji Plus Owner:

Thank you for purchasing the **Gojji Plus** Multi-Functional Monitoring System. This manual provides important information to help you to use the system properly. Before using this system, please read and follow the instructions in this Owner's Manual.

The system measures both blood glucose and  $\beta$ -ketone levels. Regular monitoring of your blood glucose and  $\beta$ -ketone levels can help you and your healthcare professional gain better control of your diabetes. Due to its compact size and easy operation, you can use the **Gojji Plus** Multi-Functional Monitoring System to easily monitor your blood glucose and  $\beta$ -ketone levels by yourself anywhere, any time.

If you have other questions regarding this product, please contact the place of purchase or call our customer service line at 1-844-GO-GOJJI (9:00 am - 6:00 pm PST, Mon. - Fri.). For assistance outside these hours, please contact your healthcare professional.

# IMPORTANT SAFETY PRECAUTIONS

## READ BEFORE USE

- The meter and lancing device are for single-patient use only. Do not share them with anyone including your family members! Do not use on multiple patients!
- All parts of the kit are considered biohazardous and can potentially transmit infectious diseases, even after you have performed cleaning and disinfection.

For more information, please visit

- “Letter to Manufacturers of Blood Glucose Monitoring Systems Listed With the FDA”  
<https://public4.pagefreezer.com/browse/FDA/12-01-2022T02:57/https://www.fda.gov/medical-devices/in-vitro-diagnostics/letter-manufacturers-blood-glucose-monitoring-systems-listed-fda>
- “CDC Clinical Reminder: Use of Fingertick Devices on More than One Person Poses Risk for Transmitting Bloodborne Pathogens” (2010)  
<https://www.cdc.gov/injectionsafety/fingertick-devicesbgm.html>
- “Infection Prevention during Blood Glucose Monitoring and Insulin Administration”  
<https://www.cdc.gov/injectionsafety/blood-glucose-monitoring.html>
- “Self-Monitoring Blood Glucose Test Systems for Over-the-Counter Use”  
<https://www.fda.gov/regulatory-information/search-fda-guidance-documents/self-monitoring-blood-glucose-test-systems-over-counter-use>

1. Use this device **ONLY** for the intended use described in this manual.
2. Do **NOT** use accessories which are not specified by the

manufacturer.

3. Do **NOT** use the device if it is not working properly or is damaged.
4. Do **NOT** use the equipment in places where aerosol sprays are being used, or where oxygen is being administered.
5. Do **NOT** under any circumstances use the device on newborns or infants.
6. This device does **NOT** serve as a cure for any symptoms or diseases. The data measured is for reference only.
7. Before using this device to test blood glucose and  $\beta$ -ketone, read all instructions thoroughly and practice the test. Carry out all the quality control checks as directed.
8. Keep the device and testing equipment away from young children. Small items such as the battery cover, batteries, test strips, lancets and vial caps are choking hazards. If swallowed, consult a doctor immediately for advice.
9. Do **NOT** use this instrument in close proximity to sources of strong electromagnetic radiation, as these may interfere with the accurate operation.
10. Proper maintenance and periodic control solution tests are essential to the longevity of your device. If you are concerned about the accuracy of your measurement, please contact local customer service or your healthcare professional for help.

**KEEP THESE INSTRUCTIONS IN A SAFE PLACE**

# TABLE OF CONTENTS

<b>BEFORE YOU BEGIN</b>	<b>1</b>
Important Information	1
Intended Use	2
Test Principle	2
Contents of System	3
Meter Overview	4
Display Screen	5
Test Strips	6
Recharging the Battery	7
Setting the Meter	7
Cellular Connectivity for the Meter	11
<b>THE THREE MEASURING MODES</b>	<b>12</b>
<b>BEFORE TESTING</b>	<b>12</b>
Calibration	12
How to Code Your Meter (for $\beta$ -ketone)	12
Checking the Code Number (for $\beta$ -ketone)	12
Control Solution Testing	13
Performing a Control Solution Test	15
<b>TESTING WITH BLOOD SAMPLE</b>	<b>17</b>
Preparing the Lancing Device for Blood Testing	17
Preparing the Puncture Site	17
Performing a Blood Glucose/ $\beta$ -ketone Test	18
<b>METER MEMORY</b>	<b>21</b>

Reviewing Test Results	21
Reviewing Blood Glucose Day Average Results	21
<b>DOWNLOADING THE RESULTS</b>	<b>23</b>
Data Transmission via Cellular Connectivity	23
<b>MAINTENANCE</b>	<b>24</b>
Battery	24
Caring for Your Meter	25
Disinfecting Procedures	26
Caring for Your Test Strips	28
Important Control Solution Information	29
<b>SYSTEM TROUBLESHOOTING</b>	<b>30</b>
Result Readings	30
Error Messages	31
Troubleshooting	33
<b>DETAILED INFORMATION</b>	<b>34</b>
Reference Values	34
<b>SPECIFICATIONS</b>	<b>36</b>

# BEFORE YOU BEGIN

## Important Information

- Severe dehydration and excessive water loss may cause readings which are lower than actual values. If you believe you are suffering from severe dehydration, consult a healthcare professional immediately.
- If your blood glucose or  $\beta$ -ketone results are lower or higher than usual, and you do not have symptoms of illness, repeat the test. If you have symptoms or continue to get results which are higher or lower than usual, follow the treatment advice of your healthcare professional.
- Use only fresh whole blood samples to test your blood glucose or  $\beta$ -ketone. Using other substances will lead to inaccurate results.
- If you are experiencing symptoms that are inconsistent with your blood glucose or  $\beta$ -ketone test results and you have followed all the instructions given in this owner's manual, contact your healthcare professional.
- The device should not be used on severely hypotensive individuals or patients in shock. Readings which are lower than actual values may occur for individuals experiencing a hyperglycemic- hyperosmolar state, with or without ketosis. Please consult the healthcare professional before use.
- **Limitation:**  
The device should not be used on individuals in hyperglycemic-hyperosmolar state, with or without ketosis; not for neonatal use; not for use on critically ill patients.

## Intended Use


The **Gojji Plus** Multi-Functional Monitoring System is intended for the quantitative measurement of glucose in fresh capillary whole blood from the finger, and for the quantitative measurement of  $\beta$ -ketone (beta-hydroxybutyrate) in fresh capillary whole blood from the finger. The **Gojji Plus** Multi-Functional Monitoring System is intended for *in vitro* diagnostic use and is intended for single-patient use as an aid to monitor the effectiveness of a diabetes control program. The system should not be used for the diagnosis of or screening for diabetes, nor for use on neonates.

## Test Principle

Your system measures the amount of sugar (glucose) or  $\beta$ -ketone in whole blood. The glucose or  $\beta$ -ketone tests is based on the measurement of electrical current generated by the reaction of glucose or  $\beta$ -ketone with the reagent of the strip. The meter measures the current, calculates the blood glucose or  $\beta$ -ketone levels, and displays the result. The strength of the current produced by the reaction depends on the amount of glucose or  $\beta$ -ketone in the blood sample.

## Contents of System

The **Gojji Plus** system kit includes:


1. Meter
2. Protective Wallet
3. Quick Start User Guide
4. Rechargeable Li-Ion Battery (installed)
5. USB Type C Cable
6. USB Adapter
7. SIM Card (installed)

Test strips, control solution, sterile lancets and lancing device may be purchased separately. Please make sure you have those items needed for a blood test beforehand.

### NOTE:

If any items are missing from your kit or opened prior to use, or to obtain test strips and control solutions, please contact the place of purchase or call Customer service for assistance.

## Meter Overview

① **Set Button (S)**

Enter setting mode or cellular connectivity mode.

② **Cellular Connectivity / Charging Indicator**

③ **Display Screen**

④ **Main Button (M)**

Enter memory mode or day average mode.

⑤ **Test Strip Indicator**

⑥ **Test Strip Slot**

Insert test strip here to turn the meter on for testing.

⑦ **USB Type C Port**

For charging the battery or data transmission.

⑧ **Test Strip Ejector**


Eject the used strip by pushing this button.

⑨ **Battery Compartment**

⑩ **Reset Button**

Reset the meter.

## Display Screen


- |                                   |                       |
|-----------------------------------|-----------------------|
| 1. Glucose Testing                | 8. Date / Time        |
| 2. Ketone Testing                 | 9. Measurement Unit   |
| 3. Test Result                    | 10. Measuring Mode    |
| 4. Error Message / Ketone Warning | 11. Volume Symbol     |
| 5. Low Battery Symbol             | 12. Test Strip Symbol |
| 6. Day Average                    | 13. Blood Drop Symbol |
| 7. Memory Symbol                  | 14. Code              |

## Test Strips

### 1. Absorbent Hole

Apply a drop of blood here. The blood will be automatically absorbed.

### 2. Confirmation Window

This is where you confirm if enough blood has been applied to the absorbent hole in the strip.

### 3. Test Strip Handle

Hold this part to insert the test strip into the meter slot.

### 4. Contact Bars

Insert this end of the test strip into the meter. Push it in firmly until it stops.

### ATTENTION:

**The front side of test strip should face up when inserting test strip.**

Test results might be wrong if the contact bar is not fully inserted into the test slot.

### NOTE:

The **Gojji Plus** Multi-Functional Monitoring System should only be used with **Gojji Plus** Multi-Functional Monitoring System Test Strips. Using other test strips with this monitoring system may produce inaccurate results.

## Recharging the Battery

Your meter comes with one 3.7V Li-Ion rechargeable battery embedded. Please recharge the meter when you get it first.


1. Connect the USB Type C cable to the meter, and connect the other end of the cable to the USB adapter.
2. Plug the usb adapter into an electrical outlet.
3. The charging indicator will light up in red, which indicates the Li-Ion battery is recharging. The meter is fully recharged when charging indicator turns green.

### NOTE:

- It may take approximately 2 hours to be fully recharged for the first time.
- Use only the USB adapter (5V/1A) provided by **Gojji**. Use of unsupported USB adapters may damage the device.

## Setting the Meter

Setting up your meter to the correct date, time, voice, etc. for better health management.

### Entering the Setting Mode

Start with the meter off (no test strip inserted). Press **S** to enter the Setting Mode.

#### 1. Setting the date

With the YEAR flashing, press **M** until the correct year appears. Press **S**.

With the MONTH flashing, press **M** until the correct month appears. Press **S**.

With the DAY flashing, press **M** until the correct day appears. Press **S**.

## 2. Setting the time format and time

Press **M** to select the desired time format (12h or 24h). Press **S**.


With the Hour flashing, press **M** until the correct hour appears. Press **S**.

With the Minute flashing, press **M** until the correct minute appears. Press **S**.

## 3. Deleting the memory

With “dEL” and “**M**” on the display, press **M** to select “no” to keep the results in memory then press **S** to skip.


To delete all the results, press **M** to select “YES”, and then press **S** to delete all the memory records.

#### 4. Setting the speaking volume

There are seven (7) speaking volume options to choose from. Press **M** until the desired speaking volume appears. To confirm your selection, press **S**.

#### 5. Choosing a language

Press **M** to select L1 / L2 / L3.

The default language for the meter is L1 English, while L2 is Spanish and L3 is Chinese.

To confirm your selection, press **S**.

#### 6. Setting the auto sent

With "AS" on the display, press **M** to select "On" or "OFF".

Press **S** to turn off the meter.

**NOTE:**

- This function is referring to the cellular connectivity. If “On” is selected, your result will be transmitted automatically right after the test.

**Congratulations! You have completed all settings!**

**NOTE:**

- These parameters can **ONLY be changed** in the settings page.
- If the meter is idle for 3 minutes during the settings mode, it will switch off automatically.

## Cellular Connectivity for the Meter

Before the first use, make sure the meter is connected to the eligible cellular service provided by supported carriers. It may take about 5 minutes to complete the cellular connectivity for the first use.

1. Press and hold **M** when the meter is off. The “PCL” then appears on the display, and the cellular connectivity indicator flashes in blue, which means the meter is connecting to the cellular signal.
2. If the cellular connectivity indicator turns solid blue, which means the connection of the meter to the cellular signal is successful.
3. The “PCL” on the display turns to “nb” or “m1”, which means the registration of the meter is complete. Press and hold **M** to turn off the meter.

### NOTE:

- Press **S** to turn off the meter if the meter registration (step 3) takes too long and is still not completed (i.e., more than 5 minutes).
- The “nb”, “m1” represent the “NB-IoT”, “LTE-M (CAT-M)” of LPWAN (Low-Power Wide Area Network) respectively, which indicates the type of cellular signal your meter connect to.

# THE THREE MEASURING MODES

The meter provides you with three modes for blood glucose testing: General, AC, and PC. You can switch between each mode by:

1. Start with the meter switched off. Insert a test strip to turn on the meter. The screen will display a flashing blood drop and “Gen”

## Voice Guide

*Please apply blood onto the test strip.*

2. Press **M** button to switch between General, AC and PC mode.

## BEFORE TESTING

### Calibration

You must calibrate the meter every time you begin to use a new vial of  $\beta$ -ketone test strips by setting the meter with the correct code. Test results may be inaccurate if the code number displayed on the monitor does not match the number printed on the strip vial.

### How to Code Your Meter (for $\beta$ -ketone)

1. Insert the code strip when the meter is off. Wait until the code number appears on the display.

#### NOTE:

- Make sure the code number on display, code strip, and test strip vial are the same. The code strip should be within the expiry date; otherwise, an error message may appear.
2. Remove the code, the display will show “OFF”. This tells you that the meter has finished coding and is ready for  $\beta$ -ketone testing.

### Checking the Code Number (for $\beta$ -ketone)

A code strip should be inserted prior using a new box of ketone test strips. Please make sure that the code number displayed on the meter

matches the number on the test strip vial before you proceed. If it matches, you can proceed with your test. If the codes do not match, please stop testing and contact Customer Service for help.

**NOTICE:**

The code number on the display is only for your reference, it may not be the actual code for your meter.

**WARNING:**

- It is important to make sure that the LCD displayed code is the same as the code on the test strip vial before testing. Failure to do so will get inaccurate results.
- If the LCD displayed code is not the same as the code on the test strip vial, and the code number cannot be updated, please contact Customer Service for assistance.

## Control Solution Testing

Our control solutions contain a known amount of glucose/ $\beta$ -ketone that reacts with test strips and that are used to ensure your meter and test strips are working together correctly.

**Do a control solution test when:**

- you first receive the meter,
- at least once a week to routinely check the meter and test strips,
- you begin using a new package of test strips,

- you suspect the meter or test strips are not working properly,
- your blood glucose/ $\beta$ -ketone test results are not consistent with how you feel, or if you think the results are not accurate,
- practicing the testing process, or
- you dropped or think you have damaged the meter.

## Performing a Control Solution Test

### 1. Insert the test strip into the meter

Wait for the meter to display the “”, “” and “**GLU**”/“**KETONE**” accordingly.

 **Voice Guide**


*Please relax during measurement.*

*Please apply blood onto the test strip.*

(For  $\beta$ -ketone Test)

*The code number is (number)*

### 2. Apply Control Solution

Shake the control solution vial thoroughly before use. Squeeze out a drop and wipe it off, then squeeze out another drop and place it on the tip of the vial cap. Hold the meter to move the absorbent hole of the test strip to touch the drop. Once the confirmation window fills completely, the meter will begin counting down. To avoid contaminating the control solution, do not directly apply control solution onto a strip.

### 3. Read and Compare the Result

After the meter counts to 0, the result of the control solution test will appear. Compare the result with the range printed on the test strip package. The result should fall within this range. If not, please read the instructions again and repeat the control solution test.

With “**QC**” displayed, the meter will store your test result in memory under “**QC**”.

## Out-of-range results

If you continue to have test results fall outside the range printed on the test strip package, the meter and strips may not be working properly. Do **NOT** test your blood, and call customer service for help.

### NOTE:

- The control solution range printed on the test strip package is for control solution use only. It is not a recommended range for your blood glucose,  $\beta$ -ketone level.
- See the **MAINTENANCE** section for important information about your control solutions.

# TESTING WITH BLOOD SAMPLE

## WARNING:

To reduce the chance of infection:

- Never share a lancet or the lancing device.
- Always use a new, sterile lancet. Lancets are for single use **ONLY**.
- Avoid getting hand lotion, oils, dirt, or debris in or on the lancets and the lancing device.
- Wash and dry your hands thoroughly after handling the meter, lancing device and test strips to prevent infection. For more information, please refer to the “**MAINTENANCE**” section.
- If the meter is being operated by a second person who is providing testing assistance to the user, the meter and lancing device should be decontaminated prior to use by the second person.

**Sharing the lancing device and lancets may increase the risk of contracting infectious diseases. Lancing device must not be used for more than one person.**

## Preparing the Lancing Device for Blood Testing

Please follow the instructions in the lancing device insert for collecting a blood sample.

## Preparing the Puncture Site

Stimulating blood perfusion by rubbing the puncture site before blood extraction significantly reduces variations between measurements.

**Please follow the suggestions below before obtaining a drop of blood:**

- Wash and dry your hands before you start.
- Select the puncture site on fingertips.
- Rub the puncture site for about 20 seconds before penetration.
- Clean the puncture site using cotton moistened with 70% alcohol and **let it air dry**.

### ► Fingertip testing

Press the lancing device firmly against the lower side of your fingertip. Press the release button to puncture your finger, then a click indicates that the puncture is complete.

#### NOTE:

- Choose a different spot each time you test. Repeated punctures at the same spot may cause soreness and calluses.
- It is recommended that you discard the first drop of blood as it might contain tissue fluid, which may affect the test result.

## Performing a Blood Glucose/ $\beta$ -ketone Test

### 1. Insert the test strip into the meter

Wait for the meter to display the “”, “” and “**GLU**”/“**KETONE**” accordingly.

#### Voice Guide

*Please relax during measurement.*

*Please apply blood onto the test strip.*

*(For  $\beta$ -ketone Test)*

*The code number is (number)*

### 2. Select the appropriate measuring mode

For glucose tests, there are three modes for measuring: Gen / AC / PC. Press **M** to select the appropriate measuring mode.

For  $\beta$ -ketone tests, there is one mode for measuring: General Test.

### 3. Obtaining a blood sample

Use the preset lancing device to puncture the desired site. Wipe off the first appeared drop of blood with a clean cotton pad.

The volume of blood sample must be at least 0.5 microliter ( $\mu\text{L}$ ) for blood glucose test, 0.8 microliter ( $\mu\text{L}$ ) for  $\beta$ -ketone test.

Gently squeeze the punctured area to obtain another drop of blood. Be careful **NOT** to smear the blood sample.

### 4. Apply the sample


While holding the meter in an angled position with the strip facing down, gently bring the strip to the sample so that the absorbent hole can sip in the blood.

Confirmation window should be completely filled if enough blood sample has been applied. Do **NOT** remove your finger until you hear a “beep” sound.

#### NOTE:

- Do not press the punctured site against the test strip or try to smear the blood.
- If you do not apply a blood sample to the test strip within 3 minutes, the meter will automatically turn off. You must remove and reinsert the test strip to start a new test.

- The confirmation window should be filled with blood before the meter begins to count down. **NEVER** try to add more blood to the test strip after the drop of blood has moved away. **Discard the used test strip and retest with a new one.**
- If you have trouble filling the confirmation window, please contact your healthcare professional or the local customer service for assistance.

#### Voice Guide

Glucose is (number)  
milligrams per deciliter.

Ketone is (number)  
millimoles per liter.

**When removing the lancet, follow the instructions from the lancing device insert.**

#### **WARNING:**

- The used lancet and test strips may be biohazardous. Please discard them properly.
- Wash and dry your hands thoroughly after handling the meter, lancing device and test strips to prevent infection. For more information, please refer to the **MAINTENANCE** section.

## 5. Read your result

The result of your test will appear after the meter counts down to 0. The result will be stored in the memory automatically.

If auto-send is on, you can press **S** after the "PCL" appears to switch off the meter, and then insert the new test strip to perform another blood test.

## 6. Eject the used test strip

Eject the test strip by pushing the eject button on the top. Use a sharps bin to dispose of used test strips. The meter will switch off automatically.

# METER MEMORY

Your meter stores the 1,000 most recent test results with date and time in its memory. To enter the meter memory, **start with the meter switched off.**

## Reviewing Test Results

### 1. Press and release **M**

“**M**” will appear on the display. The first reading you see is the last testing result along with date, time and the measurement mode.


#### Voice Guide

*Glucose is (number) milligrams per deciliter.*

*Ketone is (number) millimoles per liter.*

**2. Press **M**** to recall the test results stored in the meter each time you press. After the last test result, press **M** again and the meter will be turned off.

## Reviewing Blood Glucose Day Average Results

### 1. Press and release **M**

When “**M**” appears on the display, keep pressing **M** for 3 seconds until the flashing “DAY AVG” appears. Release **M** and then your 7-day average result measured in general mode will appear on the display.

**2. Press M to review** 14-, 21-, 28-, 60- and 90- day average results stored in each measuring mode.

### **3. Exit the meter memory**

Keep pressing **M** and the meter will turn off.

#### **NOTE:**

- Any time you wish to exit the memory, keep pressing the **M** button for 5 seconds or leave it without any action for 3 minutes. The meter will switch off automatically.
- Control solution results are **NOT** included in the day average.

# DOWNLOADING THE RESULTS

## Data Transmission via Cellular Connectivity

You can transmit your measurement data from the meter via cellular connectivity signal to view test results on your personal computer.

### 1. Setting the auto send on.

To turn the auto send on, please see the section “Setting the Meter”.

\*You can also initiate data transmission **manually** by pressing and holding the **S** button **with the meter is off**.

### 2. Data transmission

After measurement, "PCL" will appear on the meter display, indicating that the meter is in cellular connectivity mode. Results will be transmitted with date and time. The meter will automatically switch off after the data transmission is completed.

#### NOTE:

- The cellular connectivity mode is activated only when “auto send” is on.
- You can press **S** to turn off the meter during data transmission. Make sure you press **S** after "PCL" appears.

# MAINTENANCE

## Battery

Your meter comes with one 3.7V Li-Ion rechargeable battery.

### Low Battery Signal

The meter will display the message below to alert you when the meter power is getting low.

1. **The “” appears along with display message:**

The meter is functional and the result remains accurate, but it is time to recharge the battery.

2. **The “” appears with E-b and “ ” :**

The power is not enough to do a test. Please recharge the battery immediately.

The charging indicator will light up in red while recharging, and it will turn green when the battery recharge is completed.

#### WARNING:

- Recharging the battery does not affect the test results stored in the meter.
- Properly dispose of the battery according to your local environmental regulations.
- The Li-Ion battery must be recharged with the **Gojji Plus** meter by connecting USB AC adapter to an electrical outlet.

## Caring for Your Meter

To avoid the meter and test strips attracting dirt, dust or other contaminants, please wash hands thoroughly with soap and water before and after use.

### **What is Cleaning and Disinfection?**

Cleaning and disinfection are different. Cleaning is the process of removing dirt (e.g. food debris, grease, dust); disinfection is the process of killing germs (e.g. bacteria and viruses).

### **When to clean and disinfect the meter?**

Clean the meter when you see any dirt on it. You should disinfect the meter at least once a week to prevent infection.

### **How to clean and disinfect the meter?**

The meter must be cleaned prior to the disinfection. Use one disinfecting wipe to clean exposed surfaces of the meter thoroughly and remove any visible dirt or blood or any other body fluid with the wipe. Use a second wipe to disinfect the meter. **Do NOT use organic solvents to clean the meter.**

We recommend for meter cleaning and disinfection you should use the disinfecting wipe/towelette from below. The following product has been shown to be safe for use with the **Gojji Plus** Multi-Functional Monitoring System.

▶ Micro-Kill+™ (Micro-Kill Plus™) by Medline  
(EPA Reg. No. 59894-10-37549)

To obtain disinfecting wipes and other information, please contact Medline at 1-800-MEDLINE (1-800-633-5463) or visit [www.medline.com](http://www.medline.com).

## Disinfecting Procedures

- 1.** Take out one disinfecting wipe from the package and squeeze out any excess liquid in order to prevent damage to the meter.
- 2.** Wipe all meter's exterior surface display and buttons. Hold the meter with the test strip slot pointing down and wipe the area around the test slot but be careful not to allow excess liquid to get inside. Keep the meter surface wet with disinfection solution for a minimum of 2 minutes for Micro-Kill+™ wipes. Follow the instructions on the package label of disinfecting wipe. Use two or more wipes if necessary.
- 3.** Remove the wipe. Allow the meter surface to dry completely.
- 4.** Discard the used wipes and never reuse them. Wash your hands thoroughly with soap and water after handling the meter, lancing device and test strips to avoid contamination.

Improper system cleaning and disinfection may result in meter malfunction. If you have a question, please contact customer service at 1-844-GO-GOJJI for assistance.

This device has been validated to withstand 10,950 cleaning and disinfection cycles using the recommended disinfecting wipe/towelette. The tested number of cycles is estimated by 10 cleaning and disinfection cycle per meter per day for 3 years, the expected life of device. The meter should be replaced after the validated number of cleaning and disinfection cycles or the warranty period, whichever comes first.

Stop using the meter if you see any signs of deterioration. For example:

- Meter cannot be turned on
- LCD display cracks or becomes cloudy
- Buttons no longer function
- Meter outer casing cracks
- Data cannot be transmitted to pc
- Color or paint/printing on housing is abnormal
- Scratches or abrasions on meter are higher than acceptable

Please contact customer service for a replacement meter if any of the signs of deterioration are noticed.

**NOTE:**

- Do **NOT** clean and disinfect the meter while performing tests.
- If the meter is being operated by a second person, the meter and lancing device should be decontaminated prior to use by the second person.
- Do **NOT** allow cleaning and disinfecting solution to get in the test slot.
- If you do get moisture in the test strip slot, wipe it away with a corner of tissue.
- Always dry the meter thoroughly before using it.
- Do not spray the meter directly with cleaning solutions especially those containing water (i.e. soapy water), as this could cause the solution to enter the case inside and damage the electronic components or circuitry.

**Meter Storage**

- Storage conditions: -4°F to 140°F (-20°C to 60°C), below 95% relative humidity.
- Always store or transport the meter in its original storage case.
- Avoid dropping and heavy impact.
- Avoid direct sunlight and high humidity.

## Caring for Your Test Strips

- Storage conditions: 35.6°F to 89.6°F (2°C to 32°C), 10% - 85% R.H. for blood glucose test strips; 35.6°F to 86.0°F (2°C to 30°C), 10% - 85% R.H. for blood  $\beta$ -Ketone test strips. Do **NOT** freeze.
- Store your test strips in their original vial only. Do not transfer to another container.
- Store test strip packages in a cool dry place. Keep away from direct sunlight and heat.
- After removing a test strip from the vial, immediately close the vial cap tightly.
- Touch the test strip with clean and dry hands.
- Use each test strip immediately after removing it from the vial.
- Write the opening date on the vial label when you first opened it.
- Discard remaining test strips after expiration.
- Do not use test strips beyond the expiration date. This may cause inaccurate results.
- Do not bend, cut, or alter a test strip in any way.
- Keep the strip vial away from children since the cap and the test strip may be a choking hazard. If swallowed, promptly see a doctor for help.

For further information, please refer to the test strip package insert.

## Important Control Solution Information

- Use only **Goji** control solutions with your meter.
- Do not use the control solution beyond the expiration date or 3 months after first opening. Write the opening date on the control solution vial and discard the remaining solution after 3 months.
- It is recommended that the control solution test be done at room temperature 68°F to 77°F (20°C to 25°C). Make sure your control solution, meter, and test strips are at this specified temperature range before testing.
- Shake the vial before use, discard the first drop of control solution, and wipe off the dispenser tip to ensure a pure sample and an accurate result.
- Store the control solution tightly closed at temperatures between 35.6°F to 86°F (2°C to 30°C). Do **NOT** freeze.

# SYSTEM TROUBLESHOOTING

If you follow the recommended action but the problem persists, or error messages other than the ones below appear, please call your local customer service. Do not attempt to repair yourself and never try to disassemble the meter under any circumstances.

## Result Readings

### Glucose test

MESSAGE	WHAT IT MEANS
Lo	Appears when your blood glucose test result is below the lower measurement limit, which is less than 20 mg/dL (1.1 mmol/L).
Hi	Appears when your result is higher than the limit of measurement, which is higher than 600 mg/dL (33.3 mmol/L).

### $\beta$ -Ketone test

MESSAGE	WHAT IT MEANS
Lo	Appears when your $\beta$ -ketone test result is below the lower measurement limit, which is less than 0.1 mmol/L.
Hi	Appears when your $\beta$ -ketone result is higher than the limit of measurement, which is higher than 8.0 mmol/L.

## Error Messages

MESSAGE	WHAT IT MEANS	WHAT TO DO
E-b	Appears when the battery cannot provide enough power for a test.	Recharge the battery immediately.
E-2	Expired code strip.	Repeat the test with a new lot of test strip.
E-U	Appears when a used test strip is inserted.	Repeat with a new test strip.
E-E E-O E-A	Problem in operation.	Repeat the test with a new test strip. If the meter still does not work, please contact the customer service for assistance.
E-F	Appears when test strip is removed while counting down.	Repeat the test with a new test strip.
E-C	Appears when the wrong code strip is inserted or other coding errors.	Make sure to insert the right code strip for the meter. Check if the code numbers on the code strip and the strip vial label are the same.
E-t	Appears when ambient temperature is below or above system operation range.	System operation range is 46.4 °F to 113 °F (8 °C to 45 °C). Repeat the test after the meter and test strip are in the above temperature range.
E-3	Problem with the voice chip.	Please contact customer service for help.

## Cellular Connectivity Error Code List:

MESSAGE	WHAT IT MEANS	WHAT TO DO
E11 E23 E28 E31-E34 E40	Hardware failure.	Please contact customer service and send your meter back for maintenance.
E12	The network is unavailable.	Please try again at a different location. If the problem persists, please contact customer service for assistance.
E16	Error when POST data.	
E14 E24	The SIM card status or other unknown error occurred.	Please try again at a different location with better signal. If the problem persists, please contact customer service for assistance.
E27	Undefined HTTP POST error.	Please try again at a different location. If the problem persists, please contact customer service for assistance.
E35	Get signal failure.	
E53	The meter is not registered in the cloud system.	Please contact the healthcare provider for assistance.

## Troubleshooting

1. If the meter does not display a message after inserting a test strip:

POSSIBLE CAUSE	WHAT TO DO
Battery exhausted.	Recharge the battery.
Test strip inserted upside down or incompletely.	Insert the test strip with contact bars end first and facing up.
Defective meter or test strips.	Please contact customer service.

2. If the test does not start counting down after applying the sample:

POSSIBLE CAUSE	WHAT TO DO
Defective test strip.	Repeat the test with a new test strip.
Sample applied after automatic switch-off (3 minutes after last user action).	Repeat the test with a new test strip. Apply sample only when flashing “  ” appears on the display.
Defective meter.	Please contact customer service.

3. If the control solution testing result is out of range:

POSSIBLE CAUSE	WHAT TO DO
Error in performing the test.	Read instructions thoroughly and repeat the test again.
Control solution vial was poorly shaken.	Shake the control solution vigorously and repeat the test again.
Expired or contaminated control solution.	Check the expiry date of the control solution.
Control solution that is too warm or too cold.	Control solution, meter, and test strips should be at room temperature (68°F to 77°F / 20°C to 25°C) before testing.
Defective test strip.	Repeat the test with a new test strip.
Meter malfunction.	Please contact customer service.

# DETAILED INFORMATION

## Reference Values

Blood glucose plus  $\beta$ -ketone monitoring plays an important role in diabetes control. A long-term study showed that maintaining blood glucose levels close to normal can reduce the risk of diabetes complications by up to 60%<sup>\*1</sup>. The results provided by this system can help you and your healthcare professional monitor and adjust your treatment plan to gain better control of your diabetes.

Time of day	Normal plasma glucose range for people <b>without</b> diabetes (mg/dL)
Fasting <sup>*2</sup> and before meal	< 100 mg/dL (5.6 mmol/L)
2 hours after meals	< 140 mg/dL (7.8 mmol/L)

<sup>\*1</sup> American Diabetes Association. Classification and Diagnosis of Diabetes: Standards of medical care in diabetes—2022 Jan; 45(supplement 1): S17-S38.

<sup>\*2</sup> Fasting is defined as no caloric intake for at least 8 hours.

The  $\beta$ -Ketone test measures Beta-Hydroxybutyrate ( $\beta$ -OHB), the most important of the three  $\beta$ -Ketone bodies in the blood. Normally, levels of  $\beta$ -OHB are expected to be less than 0.6 mmol/L<sup>1</sup>.

$\beta$ -OHB levels may increase if a person fasts, exercises vigorously or has diabetes and becomes ill. If your  $\beta$ -Ketone result is “Lo”, repeat the  $\beta$ -Ketone test with new test strips. If the same message appears again or the result does not reflect how you feel, contact your healthcare professional. Follow your healthcare professional’s advice before you make any changes to your diabetes medication programme. If your  $\beta$ -Ketone result is between 0.6 and 1.5 mmol/L, this may indicate development of a problem that could require medical assistance. Follow your healthcare professional’s instructions. If your  $\beta$ -Ketone result is higher than 1.5 mmol/L, contact your healthcare professional promptly for advice and assistance. You may be at risk of developing diabetic ketoacidosis (DKA).

<sup>1</sup> Wiggam MI, O’Kane MJ, Harper R, et al. Treatment of diabetic ketoacidosis using

normalization of blood 3-hydroxybutyrate concentration as the endpoint of emergency management. *Diabetes Care*. 1997; 20(9): 1347-52.

**Please consult your doctor to determine a target range that works best for you.**

# SPECIFICATIONS

**Model No.:** Gojji Plus

**Dimension:** 110 (L) x 57 (W) x 25 (H) mm

**Weight:** 71 g (without battery)

**Power Source:** One 3.7V Li-Ion rechargeable battery

**Display:** LCD with backlight

**Memory:** 1,000 measurement results with date and time

Automatic detection of electrode insertion

Automatic reaction time countdown

Auto turn-off after 3 minutes without action

Temperature warning

**Operating Condition:**

50 °F to 104 °F (10 °C to 40 °C), 10% to 85% R.H. (non-condensing)

**Meter Storage/Transportation Conditions:**

-4 °F to 140 °F (-20 °C to 60 °C), below 95% R.H.

**Glucose Test Strip Storage / Transportation Conditions:**

35.6 °F to 89.6 °F (2 °C to 32 °C), 10% to 85% R.H. (non-condensing),  
up to 21 months for unopened vial.

**$\beta$ -Ketone Test Strip Storage / Transportation Conditions:**

35.6 °F to 86.0 °F (2 °C to 30 °C), 10% to 85% R.H. (non-condensing),  
up to 18 months for unopened vial.

**Glucose Measurement Units:** mg/dL

**Glucose Measurement Range:** 20 to 600 mg/dL (1.1 to 33.3 mmol/L)

**$\beta$ -Ketone Measurement Units:** mmol/L

**$\beta$ -Ketone Measurement Range:** 0.1 to 8.0 mmol/L

This device has been tested to meet the electrical and safety requirements of: IEC/EN 61010-1, IEC/EN 61010-2-101, IEC/EN 61326-1, IEC/ EN 61326-2-6, IEC/EN 60601-1-2.

For more detailed instructions, please visit [www.gojji.com](http://www.gojji.com)

## **FEDERAL COMMUNICATIONS COMMISSION (FCC) STATEMENT**

**This device complies with Part 15 of the FCC Rules. Operation is subject to the following two conditions:**

1. This device may not cause harmful interference and
2. This device must accept any interference received, including interference that may cause undesired operation of the device.

### **15.21**

You are cautioned that changes or modifications not expressly approved by the part responsible for compliance could void the user's authority to operate the equipment.

### **15.105(b)**

#### **Federal Communications Commission (FCC) Statement**

This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to part 15 of the FCC rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/TV technician for help.

#### **FCC RF Radiation Exposure Statement:**

1. This transmitter must not be co-located or operating in conjunction with

any other antenna or transmitter.

2. This equipment complies with FCC RF radiation exposure limits set forth for an uncontrolled environment. This equipment should be installed and operated with a minimum distance of 20 centimeters between the radiator radiation source and your body.

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**Distributed by Gojji, Inc. / Distribuido por Gojji, Inc.**